

## History for:

 *MMLabsCulture1\_TAGGCATG-AAGGAGTA\_L001\_R1\_001 (paired) summary report*

### De Novo Assembly

28 Nov 2015 11:05:00

Modified by: kevinmc

Add conflict annotations No

Conflict resolution Vote

Create report Yes

Create list of un-mapped reads No

Map reads to contigs Yes

Update contigs No

Minimum contig length 500

Automatic word size Yes

Match mode ignore

input:

MMLabsCulture1\_TAGGCATG-  
AAGGAGTA\_L001\_R1\_001 (paired)

Similarity 0.98

Length fraction 0.9

Insertion cost 3

Deletion cost 3

Mismatch cost 2

Color space alignment No

Global alignment No

Override paired distance No

Min distance 100

Max distance 800

Comments: 5,779,695 reads were not mapped Word size: 22

Originates from:



MMLabsCulture1\_TAGGCATG- AAGGAGTA\_ L001\_R1\_001 (paired)